c-prodirect

# Non-Ambulatory Abduction Dorsiflexion Mechanism Handbook

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# English

## FITTING INSTRUCTIONS

**\*Please read in conjunction with fitting instruction insert.**

### **Note:**

The ADM is fitted with the spring mechanisms located on the outside of the leg and the sandal buckles located on the inside of the foot.

- 1. Place the foot into the open sandal with the heel back.**
- 2. Keeping the heel back, fasten the top buckle strap firmly.**
- 3. Fasten the top strap around the ankle.**
- 4. Check the heel is still back and down by looking through the holes in the back of the sandal. If not, go back to instruction 2.**
- 5. Fasten the other strap(s) ensuring toes are straight and flat in the sandal. These do not need to be tight.**
- 6. Clip the ADM into the sandal.**
- 7. Place the ADM as far around the leg as it will go. This engages the TTJ spring and aligns it into the correct position.**
- 8. Whilst holding the ADM in place, fasten the ADM straps. The lower strap should be fastened first and firmly to keep the ADM in position. The upper strap should be fastened secondly and does not need to be as firm as the lower strap.**
- 9. Press down on the straps to secure them.**
- 10. Cut any excess material using the guides shown.**
- 11. Once the child is settled and relaxed in the ADM abduction and dorsiflexion should be visible. If you are unsure that the correct level of abduction or dorsiflexion is being achieved, consult your clinician or see our troubleshooting guide.**
- 12. The ADM comes preset. Do not adjust the ADM or tamper with it in any way.**

## INTRODUCTION

The Abduction Dorsiflexion Mechanism (ADM) is a dynamic Ankle Foot Orthosis that assists anatomically correct foot dorsiflexion, foot abduction and foot eversion motions.


What is the ADM and what does it do?

### Ankle Foot Motions


The two major joints of the foot are the Tibio-Talar Joint (TTJ) and the Sub-Talar Joint (STJ).

TTJ: The dominant motion is Dorsiflexion / Plantarflexion

8 degrees from Transverse plane

20 to 30 degrees from frontal plane

STJ: The dominant motions are abduction /adduction and eversion/inversion

42 degrees from transverse plane


16 degrees from sagittal plane

### ADM Design

The complete ADM unit is based on the two joints above and comprises of a sandal or shoe and the Abduction Dorsiflexion Mechanism (ADM).

The ADM TTJ (dorsiflexion) & STJ (abduction) mechanism axes are powered by torsion springs.

The ADM TTJ mechanism axis is aligned to the mean axis of motion of the Tibio-Talar Joint.

The ADM STJ mechanism axis is aligned to the mean axis of motion of the Sub-Talar Joint.

### Indications

The Non-Ambulatory ADM is intended for use on successfully corrected clubfeet in accordance with the Ponseti method after completion of casting. It may also be used at later stages of treatment. It can also be used for other conditions which affect dorsiflexion and abduction of the foot such as neurological conditions. Use of the ADM should be as directed by a clinician.

### How the ADM works

The Abduction Dorsiflexion Mechanism (ADM) features two spring mechanisms. Each ADM (mechanism) contains a torsion spring with the lower spring mechanism (STJ) responsible for abduction and the upper spring mechanism (TTJ) responsible for dorsiflexion.

The ADM stretches tendons and ligaments and allows the child to move their foot to maintain flexibility and condition. These stretches will be most efficient when the child is in a relaxed state or sleeping.

### Intended Use

This should always be as directed by a clinician.

In usual circumstances for clubfoot, the brace can be used for the standard protocol of 23 hours a day for 3 months, followed by 12-14 hours of nap or night time use after completion of casting in accordance with the Ponseti method.

For other conditions, protocol may vary.

### Choosing the Spring Strength

The spring strength should always be determined by a clinician.

The guidelines for use are:

For clubfoot, standard spring options are suitable for most patients with well corrected feet.

If standard spring strengths do not provide the required abduction and dorsiflexion a clinician may advise softer or stronger springs.

### Introducing the child to an ADM

There may be an adjustment period when first starting with the ADM. If the child is struggling with the time recommended by your clinician, we suggest that you gradually build up the time over a period of 1-2 weeks.

### TROUBLESHOOTING

We have included some common queries below to help you troubleshoot with the Non-Ambulatory ADM.

- How do I know which spring strength I have? For each size sandal, the corresponding ADM will have specific spring strengths. To find the spring strength of the ADM you are using, pull back the material at the top of the lower leg attachment to show the spring strength label. The label should look similar to this: **STJ 90 TTJ 120**  
This specific label is for an abduction spring strength (STJ) of 90 and a dorsiflexion spring strength (TTJ) of 120.
- How do I know which sandal size I have? The size of the sandal is printed underneath the sandal on the black sole.
- How do I know which ADM height I have? To check which ADM height you have received, look at the back of the lower leg attachment under the logo. For example "ADM SMALL" is a size small ADM. The ADM is available in six sizes: Extra Small, Small, Small-Medium, Medium, Medium-Large and Large dependent on the sandal size it corresponds to.
- Should the child wear socks? Make sure the child is wearing close fitting, smooth socks that extend above the ADM. Badly fitting socks may cause sores.
- Can the child walk and play in the Non-Ambulatory ADM? The Non-Ambulatory ADM is not designed for excessive

physical activity, so we recommend children only wear the ADM indoors and for occasional walking in the house, such as going to bed and going to the bathroom.

- **The child is getting red marks or sores:** This could be happening for a number of reasons:
  - Incorrect size sandal and/or ADM - please check with a clinician.
  - Incorrect fitting of the sandal and/or ADM - please refer to our Fitting Instructions.
  - Padding interfering with foot motion - please remove the padding.

Sensitive skin may still get sores or redness underneath the sandal straps. It may help in these circumstances to replace the tongue provided with an alternative form of padding such as neoprene or sheepskin.
- **Should we use Pressure Saddles?** We do not recommend using pressure saddles or any other type of padding except as indicated above. The sandal straps should not be overlapping. In our experience the friction of extra padding may cause sores. We have also found adding these may affect the spring motions, stopping them from moving to their fullest range. If the child is getting sores, see above.
- **The ADM isn't moving the foot:** You may have found that when you fit the mechanism that the foot is not abducting or dorsiflexing as much as it should. Reasons for this can include:
  - The foot is resisting the ADM. It is absolutely vital to check the degree of abduction and dorsiflexion when the child is fully relaxed. The ADM is a dynamic brace that allows the foot position to change.
  - Wrong size sandal and corresponding ADM
  - ADM and/or sandal not fitted correctly
  - Different ADM spring strength is required
  - ADM has slipped round the leg
  - Not all feet are as flexible as others. The ADM should take the foot to where it can naturally go.

If you are concerned please speak to a clinician.
- **The ADM is squeaky:** If your ADM is squeaking, you can use a natural oil such as coconut oil to lubricate the mechanism. If it continues, please contact your supplier for a replacement.
- **What readings should I expect to see on the ADM dials once it is fitted properly and my child is relaxed?** Most doctors agree 10+ degrees on TTJ (dorsiflexion) and 15+ degrees on STJ (abduction), although it does depend on the child and the condition being treated. Some feet are more or less stiff than others particularly if they are hyper mobile, atypical, complex, or have other conditions affecting mobility etc. The ADM will only be able to move the foot within its available range of motion. Please speak to a clinician if you have any concerns.
- **When do I need a new ADM and sandal?** The ADM and sandal should be changed when they begin to become less effective i.e. when less dorsiflexion and/or abduction is achieved than when first fitted. This may correspond to the toes reaching the edge of the sandal but may be sooner. It is advisable to take a photo of the ADM when it is first fitted and providing optimal abduction and dorsiflexion and use this as a reference point to ensure the ADM is continuing to be effective. Each time the sandal size increases the mechanism will also need to be changed to accommodate leg growth and the extra force required to move the foot.
- **I have bought a bigger ADM and sandal to accommodate growth:** You should not buy an ADM and sandal bigger than required as this may not be comfortable for the child. The foot may slip in a large sandal, which may not maintain correction, and the ADM springs may be too strong for the foot and cause the child discomfort and sores. If you need a smaller size, please contact your supplier for an exchange.

## DIFFERENCES BETWEEN THE ADM AND THE BOOTS AND BAR

	Boots and Bar	ADM
Adjustments	Angle of external rotation on each side, length of bar	<b>NONE REQUIRED:</b> Comes preset with springs appropriate for patient (can be adjusted by manufacturer only)
Change to next size	When toes reach end of sandal	When toes reach end of sandal OR when ADM does not provide appropriate abduction and dorsiflexion
Padding	Suitable for pressure saddle and other padding	Padding can get in the way of the foot dorsiflexing. We recommend thicker socks instead or replace tongue with sheep skin or similar
Socks	Yes	Yes
Non clubfoot braced	Yes	No
Minimum foot length	6cm	8.1cm (smaller available by exception)
Minimum length from knee crease to bottom of heel	N/A	10.5cm
Maximum foot length	20.2cm	Unlimited but sandal/shoe may vary. Available for older children and teenagers
Growing Room	Provided in sandal, can go up a size so long as foot held in good position in sandal	Provided in sandal, do not go up a sandal size as the corresponding ADM springs may then be too strong

## STORAGE AND CLEANING

Please store the ADM and sandal out of direct sunlight or heat. Do not put the ADM in the washing machine. The ADM device should not be submersed in water but can be wiped clean with a damp cloth and dried. The ADM sandal can be washed in a washing machine on a cool gentle cycle using mild detergent. It should be left to dry naturally and a suede (or nail) brush can then be used on the straps to soften them.

## FOR CLINICIANS

If you are a medical professional, please note we do offer training to help you with the fitting of the ADM. If this is something you would like, please contact our Clinical Support team at [clinicalsupport@c-prodirect.co.uk](mailto:clinicalsupport@c-prodirect.co.uk) or your ADM supplier for more information.

### Fitting Stocks

We also offer medical professionals a free Fitting Stock to help choose the correct spring strength. This is made up of a sample stock of each size sandal and left and right ADMs of each standard spring strength. You can select the correct sandal size and the correct spring to order in the comfort of your clinic. For more information, please contact the Clinical Support team at [clinicalsupport@c-prodirect.co.uk](mailto:clinicalsupport@c-prodirect.co.uk) or your ADM supplier.

## WARNINGS

Use only in accordance with these instructions and as directed by a clinician.

Failure to use ADMs with sufficient spring strength for the patient may result in a loss of correction.

Do not dis-assemble or tamper with the ADM or ADM Sandal. All product warranties are voided for ADMs that have been tampered with in any way.

Do not use the ADM if any part is damaged, not functioning or does not correctly fit the patient.

## WARRANTY

We aim to provide the highest service levels for all of our customers. If the merchandise you received is defective or not as you ordered, please contact us for a quick resolution.

- Broken ADM (within 12 months of purchase): A free of charge replacement will be provided if an ADM has a manufacturing defect or is damaged through normal use within 12 months of purchase.
- Broken ADM Sandal (within 6 months of wear): A free of charge replacement will be provided if an ADM Sandal has a manufacturing defect or is damaged through normal use within 12 months of purchase.
- Broken ADM Sandal (after 6 months of wear): If a sandal breaks or becomes damaged through normal use and has been in use for more than 6 months then we require a new foot measurement to be taken to ensure that the current size is still appropriate for the child. If the current size is still appropriate a free of charge replacement will be provided.
- Incorrect size or incorrect ADM Springs: Please notify us within 7 days of receipt. Provided the product is returned in a saleable condition then a replacement of the correct size or configuration will be provided free of charge.
- Unwanted Product: If you are dissatisfied a refund will be provided if we are notified within 7 days of receipt and the product is returned in a new, saleable, unmarked condition and in original packaging. The products must be returned within 30 days of invoice date.

## RESOURCES

### Website

Please see our website for resources such as Downloads, Videos, Case Histories and Customer feedback:

<https://www.c-prodirect.co.uk/resources>

### ADM Facebook Page

There is a Facebook page called ADM Brace - a forum to provide questions, answers and experiences with the Abduction Dorsiflexion Mechanism (ADM). If you would like to join, follow this link and answer the questions to be approved:

<https://www.facebook.com/groups/374379346246818/>

### Clinical Support

Regarding ADM ordering, sizing, fitting and troubleshooting you can contact our Clinical Support team at [clinicalsupport@c-prodirect.co.uk](mailto:clinicalsupport@c-prodirect.co.uk). C-Pro Direct does not under any circumstances provide medical advice.

### YouTube

Our latest videos are posted on our YouTube channel. You can follow the channel here:

[https://www.youtube.com/channel/UC7YP5-UInsR5IBn\\_4uDbzRw](https://www.youtube.com/channel/UC7YP5-UInsR5IBn_4uDbzRw)

## PATENTS

The ADM is developed and manufactured by C-Pro Direct Ltd in the UK. The ADM is Patented in many jurisdictions worldwide and Patent applied for in the USA. For Europe EP Patent No EP2637612 and Registered Community Design No 002238881-0001 apply.

## GLOSSARY

- Non-Ambulatory: Unsuitable for walking in (see Troubleshooting).
- Foot Abduction: The movement of the foot away from the midline of the body.
- Foot Dorsiflexion: The movement of the foot in which the foot flexes in an upward direction.
- Foot Eversion: The combination of the foot dorsiflexing and abducting.
- Plantar Flexion: The movement of the foot in which the foot flexes downward toward the sole.
- Ponseti Method: A method of treating clubfoot. For more information, see [www.ponseti.info](http://www.ponseti.info).

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