

ADM Size Details (March 2020)

L Foot Length Range (mm)	ADM Size	Sandal Size	H2 Max Height (mm)	H1 Min Height (mm)	Lower Strap Circumference Range (mm)	Upper Strap Circumference Range (mm)
81 to 92	XXXS	0 to 1	105	55	100 to 200	120 to 220
93 to 102	XXS	2 to 3	120	56	80 to 180	100 to 200
103 to 122	XS	4	148	76	80 to 180	110 to 210
123 to 132	S	5	176	81	80 to 180	130 to 230
133 to 142	SM	6	196	93	80 to 180	130 to 230
143 to 162	M	7 to 8	220	106	90 to 200	140 to 250
163 to 172	ML	9 to 10	223	103	90 to 210	150 to 270
173 to 202	L	11 to 12	248	103	100 to 240	150 to 290

Notes:

- The top of the Leg Attachment (**H2**) should be between the upper extents of the main calf muscle and just below the popliteal crease
- The lower end of the Leg Attachment (**H1**) should be above the ankle joint
- The straps fit most children and the straps may be trimmed to remove excess length. If necessary custom length straps are available. Please contact us if you require assistance
- Heights **H1** and **H2** are based on the ADM being fitted in the neutral condition as recommended

