

# Dobbs ADM<sup>®</sup> Fixed External Rotation Bar System Overview

c-prodirect

An ultra low weight, comfortable and easy to manage fixed external rotation brace system for the prevention of clubfoot relapse.



C-Pro Direct Ltd, Edenbridge, Kent, TN8 6HF United Kingdom.

C-Pro Direct Ireland Limited, 77 Camden Street Lower, Dublin, D02 XE80. Ireland

## The Dobbs ADM® Fixed External Rotation Bar System

Used with the ADM® Ankle Foot Orthosis (AFO) the Dobbs ADM® Fixed External Rotation Bar system forms a traditional clubfoot brace (Boots and Bar) for the prevention of clubfoot relapse.

### Indications:

- Maintenance of correction of well corrected clubfeet as per the Ponseti Method
- Maintenance of vertical talus correction
- Correction of metatarsus adductus

### Advantages and Benefits:

- Easily and rapidly configured to the clinical requirement.
- Zero abduction with zero dorsiflexion clips available for Vertical Talus.
- No screws or adjustable parts reduces complexity for clinical staff, parents and carers and assures the required clinical configuration is maintained.
- Ultra low weight (typically over 40% lighter) compared with popular alternatives.
- Attractive modern sleek appearance appearance.
- Cost and time efficient for clubfoot clinics.



Select the required angle Left and Right external rotation bar clips and the required length bar. Then assemble brace in just a few seconds.



*The ADM® Boots and Bar system is so easy to fit, intuitive to use and quick to configure. We now spend half the time in clinic fitting clubfoot braces and training parents and carers. Quote from a leading UK clubfoot clinic lead*